#rethinking critical care

hygh-tec® odor-solution

Stuhldrainage-System | Stool drainage system |
Système de drainage fécal | Sistema di drenaggio delle feci | Faecaal drainagesysteem |
Ulosteenhallintajärjestelmä

hygh-tec® bag

Stuhlsammelbeutel | Stool collection bag |
Poche de recueil de matières fécales | Sacca di raccolta delle feci | Ontlastingsopvangzak |
Ulosteenkeräyspussi

Gebrauchsanweisung | Instructions for use |
Mode d'emploi | Istruzioni per l'uso | Gebruiksaanwijzing | Käyttöohje

Deutsch | English | Français | Italiano | Nederlands | Suomi

hyghtec.com

hygh-tec®

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Dear patient, dear user,

According to the European Regulation on Medical Devices (2017/745), Advanced Medical Balloons GmbH is obliged to make you aware of the following: All serious incidents that occur in connection with the product must be reported to Advanced Medical Balloons GmbH and the competent authority of the member state in which the user and/or the patient is located.

Note:

Serious incident means the death of a patient, user or other person or the temporary or permanent serious worsening of the health of a patient, user or other person. It does not matter whether this occurred or could occur.

You can find the contact details for the competent authority of your member state online using the search terms "Competent Authorities for Medical Devices EU".

Contact details for Advanced Medical Balloons GmbH

quality@amb-medtec.com, +49 7254 40397-0

Enquiries about product training

training@amb-medtec.com, +49 7254 40397-0

Content

Product description and clinical benefit	18
structural design	19
function	20
purpose of use, indications and contraindications	20
risks of use	22
precautions	24
use of the product	25
product content	30
Overview of products	86

hygh-tec® odor-solution | Product description and clinical benefit

hygh-tec® odor-solution is the odour-reducing variant of hygh-tec® basic-plus. Unlike silicone-based stool drainage systems, hygh-tec® odor-solution has a drainage head made of polyurethane that positions the drainage system past the sphincter muscle and seals it. The balloon on the drainage head is dumbbell-shaped and has already been moulded to its full working dimensions during production.

The trans-anal balloon segment [B] is positioned inside the anal canal. An elastically deformable, trans-anal shaft element [C] runs through the inside of the drainage head. This shaft element takes up the stool in the rectum and drains it through the anus. The stool-draining trans-anal shaft element [C] has a wave-shaped profile, which supports its elastic straightening. When the sphincter muscle tone is normal, the shaft element folds into the anal canal in a radial direction. If the sphincter muscle tightens, the shaft element in the anal canal straightens in an elastic manner, releasing the drainage pathways for the stool to be drained.

Unlike silicone-based stool drainage systems, the drainage balloon is inflated with air, not water. The quantity of air to be used is specified and makes the balloon [A + B + D] tension-free and flaccid. The loosely filled balloon envelope clings in a sealing manner to the respective anatomy of the rectum and anus. The anterior, intra-rectal segment of the balloon [A] absorbs the respective force in the rectum and at the same time translates it into a sealing of the anal canal in the central, trans-anal segment of the balloon [B]. The trans-anal segment of the balloon [B] therefore adapts to the respective tone and the respective level of opening of the sphincter muscle. This enables continuous sealing of the anus even in patients who are conscious and mobilised.

The drainage head is inflated with an inflating syringe included with the product [N]. In order to do this, the syringe is connected to the inflating port [J] labelled "Air". The volume to be used for the inflation (85 ml of air) is marked on the cylinder of the inflation syringe.

Another port [K] marked "Irrig." is used to supply liquid to irrigate the system to clean it or for the rectal irrigation of the patient.


The drainage tube [G] connects to the drainage head [A – E] positioned transanally and connects the head unit [A – E] to the stool collection bag [M].

The connection to the bag is by means of a bayonet fastening. The extraction port for stool sampling [I] is found in the lower third of the drainage tube. The drainage tube [G] can be connected or blocked at any height using a freely adjustable closure strap [L].


The stool collection bag [M] is equipped with a carbon filter to release intestinal gas. The contents of the bag is hidden from the user by means of a large-scale print on the bag. The stool collected in the bag can be viewed and inspected by the user on the transparent reverse of the bag.

hygh-tec® odor-solution | structural design

A – intra-rectal balloon segment | B – trans-anal balloon segment |
 C – trans-anal shaft element | D – pre-anal balloon segment |
 E – gel-olive | F – yellow ring for position control | G – drainage tube |
 H – irrigation line | I – sampling port | J – inflation port of the drainage head |
 K – irrigation port | L – closure strap | M – stool collection bag with connector |
 N – inflating syringe

hygh-tec® odor-solution | function

hygh-tec® odor-solution dynamically adapts to the respective position of the anal sphincter muscle. The loosely inflated drainage head takes up the pressure in the rectum in the intra-rectal balloon segment (1) and uses it to seal the sphincter muscle in the trans-anal balloon segment (2) at the same time.

The trans-anal shaft element (3) by hygh-tec® odor-solution folds up elastically in the anal canal when the sphincter muscle done is normal and straightens in the respective opening of the anus in an elastic manner if the tone drops. hygh-tec® odor-solution therefore helps with the free and unhindered flow of stool.

hygh-tec® odor-solution | purpose of use, indications and contraindications

hygh-tec® odor-solution is the odour-reducing variant of hygh-tec® basic-plus. hygh-tec® odor-solution is used for the continuous, transanal drainage of stools that will flow and liquid stools.

The hygh-tec® bag stool collection bag is used to collect stool after drainage through the catheter within a stool drainage system.

Both products are used in professional healthcare facilities, for example in intensive care. The products are used by intensive care staff or professional, trained nursing staff.

The products are for patients over the age of 18 who are critically ill and in intensive care, immobile or immobilised patients and/or in some cases mobile patients with persistent diarrhoea.

hygh-tec® odor-solution and hygh-tec® bag are intended for single use.

The following applies when using the product hygh-tec® odor-solution:

- continuous, transanal drainage of stools that will flow and liquid stools over a period of not more than 30 days.

The indication for the use of the product is made exclusively by the doctor treating the patient. The type and severity of the underlying diseases and the risks associated with the use of the product are weighed against the expected benefit of continuous stool drainage by the doctor when making the indication.

The following are contraindications for the use of the product hygh-tec® odor-solution:

- Unclear, not diagnostically clarified, traumatic, ischaemic or inflammatory damage to the rectum or the anus.
- Unclear, not diagnostically clarified mass in the region of the rectum or the anus.
- Unclear, not diagnostically clarified bleeding in the region of the rectum or the anus and accumulations of blood in the patient's stool.
- Diagnostically clarified, malignant mass in the region of the rectum or the anus.
- Previous surgery in the patient's colon, rectum or anus in the last six months.
- Phase of spinal shock.
- Known allergic hypersensitivity to the materials used in the product (polyurethane, PVC, polycarbonate, silicone).
- Use in patients under the age of 18.
- Use in pregnant patients.
- Formed stools or stools that will not flow.
- Use of the product over a period of more than 30 days.

There are no indications, contraindications and restrictions in patient population specific to the hygh-tec® bag stool collection bag. Information about the indication, contraindication and patient population are based on the specifications for the stool drainage system.

hygh-tec® odor-solution | risks of use

Warnings

The following warnings must be taken into account when making the indication for the use of the product and over the course of the use of the product:

- Known, diagnostically clarified, chronically degenerative, inflammatory or ischaemic changes to the rectum and the anus must be assessed by the doctor providing treatment and ruled out as a potential contraindication for the use of the product.
- Previous surgery in the patient's colon, rectum or anus in the last six months must be clarified by the doctor providing treatment in terms of the postoperative conditions and postoperative perfusion situation before the product is used.
- Strictures or stenoses in the region of the rectum or the anus that present a potential obstacle when inserting and positioning the drainage must be medically assessed before the product is used.
- Advanced haemorrhoids that tend towards bleeding must be medically assessed before the product is used.
- In the case of patients with a known tendency to bleed and those on anticoagulant therapy, it is necessary to watch for potential accumulation of blood in the drained stools over the course of use of the product.
- The respective risk of lesions of the perianal skin and the skin on the buttocks and the upper thigh caused by pressure must be assessed in patients with disorders of rectal, anal or pre-anal sensitivity. This applies in particular in the case of patients with injuries to the spinal cord.
- Treatment with catecholamines to support the circulatory system: catecholamines are regularly used to stabilise the circulation in patients who are critically ill. Catecholamines can reduce circulation (perfusion) in the peripheral tissues of the body to a critical extent and therefore result in the development of tissue damage caused by perfusion. Damage of this type can occur in the rectum in particular during continuous stool drainage, with the intra-rectal segment of the balloon lying on the mucous membrane in the rectum over extended periods of time. It is necessary to check for potential ulcerations of the rectal mucous membrane over the course of stool drainage, particularly in the case of the administration of high doses of catecholamines with a vascular narrowing, vasoconstrictive effect. Higher grade damage to the rectum can be identified by the user by potential accumulation of blood in the stool.

- The balloon on the drainage head may only be inflated with air.
- Mixing up the inflation port on the drainage head with the irrigation port: Both ports have a Luer connector. The correct connection to the irrigation port of the drainage must therefore be ensured, particularly when irrigating the system to clean it using a syringe filled with water. The port for the supply of liquids is marked with the letters "Irrig". The port for inflating the drainage head with air is marked "Air". Accidental injection of liquids into the drainage head leads to an increase in pressure in the drainage balloon. This can result in the tube element of the drainage head which drains the stool collapsing, permanently closing the drainage and resulting in subsequent abdominal congestion of stools and intestinal gas.
- If stools are discharged from or leak out of the patient's anus over the course of the use of the product, the balloon in the drainage head should be deflated completely and then re-inflated with 85 ml of air. Inflating above 85 ml does not improve the seal. Overinflating the balloon can lead to a partial or complete closure of the tube components draining the stools inside the drainage head. Permanently closing the lumen draining the stool can lead to the abdominal congestion of stools and intestinal gas.
- If the balloon of the drainage head is deflated and re-inflated with 85 ml of air, it is necessary to ensure that the balloon is fully deflated before re-inflating. This means overinflation of the balloon caused by the residual volume remaining can be avoided. Overinflating the balloon can lead to a partial or full closure of the lumen of the drainage which drains the stools.
- When infusing or injecting liquids via the irrigation port on the drainage, no substances that are harmful to the mucous membrane may be applied. The indication for the rectal application of substances may only be made by the doctor providing treatment. Any toxic interactions between the substance administered and the drainage materials (polyurethane, PVC, polycarbonate, silicone) must be clarified before application.
- No media with too high a temperature that harm the mucous membrane may be supplied when irrigating the system to clean it and during rectal irrigation.
- Formed, partially formed or thick, mushy stools cannot be drained via the drainage. A continuous flow of stools must therefore be ensured during use and the patient checked for the abdominal congestion of stools and intestinal gas. If the stools are insufficiently able to flow, use of the product must be discontinued.

hygh-tec® odor-solution | precautions

- The patient's anal canal and rectum must be probed or digitally inspected before insertion of the drainage. Any masses or signs of rectal or anal bleeding must be diagnostically clarified. Formed, partially formed or thick, mushy stool residue in the rectum should be removed where possible by means of the rectal insertion of a conventional irrigation system before the product is inserted.
- When digitally probing the rectum, the tone of the anal sphincter muscle should be assessed. If there is no or insufficient tone or the anus is gaping open, the efficacy of the seal and the anchoring of the drainage may be impaired. If the tone of the sphincter muscle is pathologically high, this must be assessed by the doctor providing treatment with regard to the likelihood of the development of damage to the anus caused by pressure.
- If the drainage is expelled from the rectum, the drainage head can be reinserted by the user. In order to do this, the drainage head is completely removed and once again inflated with 85 ml of air following further transanal positioning.
- To prevent ischaemic damage to the rectum caused by pressure, potentially resulting in ulceration or necrosis, it is necessary to ensure that the specified filling volume of 85 ml of air for the drainage head is complied with and the stool drainage tube is not placed under persistent tensile stress for the entire duration of use.
- Inform the doctor providing treatment immediately if the patient indicates any pain or irritation in their abdomen, rectum or anus.
- The stools drained must be checked for potential accumulation of blood over the course of the use of the stool drainage. If this is observed, potential causes related to a lack of perfusion must be clarified by the doctor providing treatment immediately.
- After long periods of drainage it can take up to two days for full sphincter muscle function to be restored.
- In the case of the infusion or injection of liquids into the patient's rectum, media that are at body temperature should be avoided to avoid unwanted rectal stool emptying reflexes.
- When using the product, particularly when irrigating the system to clean it, it is necessary to ensure a sufficient gradient between the patient and the collection bag placed below the patient.
- Possible changes in the patient's fluid balance must be taken into account when supplying liquids rectally in patients with heart failure and kidney failure.
- The drainage system must not be used in conjunction with other catheters or devices positioned in the anal canal.

- The patient's body should ideally not be laid on the drainage tube to prevent damage to tissue caused by pressure.
- Persistent discharge or leakage of stool in the region of the patient's perianal and sacral skin can cause irritation and inflammatory reactions of the epidermis there.
- The product may only be used if the packaging is not damaged.
- The product is designed for single use and only intended to be used on a single patient. The products and the components attached to it may not be processed and used in another patient.
- In the case of obese patients, it is necessary to check for damage to the skin in the buttock and upper thigh regions caused by pressure. This can be caused by the pre-anal part of the drainage head that is positioned there.
- The continuous drainage of stools with a particularly high fat content can lead to the partial stiffening of the drainage tube that provides drainage. In these cases, the drainage system should be removed and replaced with a new drainage system.

hygh-tec® odor-solution | use of the product

A – Preparing the patient

1. Position the patient lying down on their left side, if possible, to insert the drainage system.
2. Palpate the rectal ampulla. Use a lubricant.
Check for potential masses or narrowing within the rectum and the anus.
3. If there are stools in the rectal ampulla, the patient should be prepared by administering a rectal enema. Use a conventional rectal irrigation set for this purpose.
4. Ensure that there is sufficient tone in the anal sphincter muscle.

B – Preparing the system

1. Connect the drainage system to the stool collection bag provided.
Twist the locking ring until it snaps into place.
2. Completely deflate the drainage system balloon using the inflation syringe provided.
3. Coat the front third of the drainage head system with lubricant gel.

C – Inserting the system

1. Insert the drainage head carefully around two-thirds of the way into the anus.
2. Inflate the drainage head with 85 ml of air. Use the inflating syringe provided. Ensure that the speed of inflation is moderate and not too quick to avoid triggering an anal

opening reflex or an elimination reflex. The balloon will be a bit underinflated and soft in the rectum. Do not exceed the specified inflation level of 85 ml of air.

3. Check the transanal position of the drainage head. The yellow position control ring must be visible outside of the anus. With normal anatomy, the lower section of the balloon part of the drainage head will protrude from the anus.
4. Check the degree of inflation of the section of the balloon protruding from the anus. The balloon must look and feel inflated.
5. Uncoil the drainage tube. Ensure that the drainage system outlet is clear and unobstructed. Attach the collection bag low enough below the patient.

We also recommend:

6. Carry out the first irrigation of the system approximately 30 minutes after inserting the drainage system. Preferably use 500 ml of a crystalloid infusion solution that is at body temperature. Use a conventional infusion set. It is important to connect the irrigation solution to the irrigation port marked "Irrig" and to ensure a moderate flow of the solution. If the flow is too quick it can trigger an anal opening reflex or an elimination reflex.

D – Drainage system maintenance – once per shift

1. Check the correct transanal position of the drainage head. The yellow ring must be visible outside of the anus.
2. Check the inflation of the drainage head. The balloon part of the drainage head protruding from the anus must be look and feel inflated. If the inflation level of the balloon is not clear or the inflation level cannot be assessed by the user, the drainage head should be fully deflated and then re-inflated with 85 ml of air. Use the inflating syringe provided.
3. Irrigate the system to clean it once per shift with 150 ml irrigation solution. Preferably use a crystalloid infusion solution and a conventional infusion set. Ensure a moderate flow of the solution. If the flow is too quick it can trigger an anal opening reflex or an elimination reflex. It is important to connect the irrigation solution to the irrigation port labelled "Irrig".
4. Ensure that the consistency of the stools is sufficiently able to flow.
5. Check the free, unimpeded flow of stool from the patient to the collection bag. The tube should be lying freely and not twisted. Monitor for signs of stool congestion. Check the patient's bowel activity on a regular basis by auscultation.

Every three days:

Completely deflate the drainage head and re-inflate the balloon with 85 ml of air. Use the inflating syringe included.

E – Additional information and instructions:

If stools are discharged from the anus:

If stools still leak from the anus despite the fact that the drainage system is correctly transanally positioned, the following causes are possible:

- rectal congestion of stool caused by a shift in the drainage lumen, in particular as a result of stools that do not flow sufficiently or cannot drain;
- overinflation of the drainage head with a full closure or collapse of the drainage lumen and subsequent rectal congestion;
- insufficient tone in the anal sphincter muscle;
- an agitated and especially vigilant patient;
- insufficient inflation of the drainage head;

If the leakage of stools is identified, the drainage head should be fully deflated and then re-inflated with 85 ml of air. Inflating above 85 ml does not improve the seal on the drainage system. Overinflating the drainage head can lead to a collapse of the tube components draining the stools inside the drainage head that closes the lumen. A closure of this type can lead to rectal congestion of stool and therefore to the system overflowing.

If the inflation level of the balloon is not clear or not able to be assessed:

The balloon on the head unit should then be fully deflated and then re-inflated with 85 ml of air. Use the inflating syringe provided.

If the balloon cannot be deflated using the inflating syringe, damage to the inflation port of the drainage head [J] may have occurred during the course of use. If this is the case, cut the inflation port off. The air will release automatically. The stool drainage system can then no longer be used. A new drainage system must be inserted to continue use.

If the drainage system slides out of the rectum:

There could be various causes of this, such as:

- rectal congestion of stool caused by a shift in the drainage lumen, in particular as a result of stools that do not flow sufficiently or cannot drain;
- insufficient tone in the anal sphincter muscle;
- an agitated and especially vigilant patient;
- insufficient inflation of the drainage head;

The drainage system can generally continue to be used. Wipe the drainage head with a dry cloth and reinsert it into the rectum according to section C "Inserting the system".

If the drainage head slides into the rectum:

The sliding of the drainage head positioned transanally in the rectum has been observed particularly in patients in a sitting position. The correct positioning of the drainage system must therefore be carefully checked in patients who are sitting and those with their upper body elevated. The correct position can be ensured by the yellow ring, which must be visible outside of the anus.

F – Rectal irrigation

1. Position the patient lying down on their left side, if possible.
2. Position the red closure strap approx. 10 cm from the patient's anus and tighten the strap.
3. Preferably use a crystalloid infusion solution at body temperature and a conventional infusion set for the rectal irrigation.
WARNING: Connect the irrigation solution to the irrigation port labelled "Irrig".
4. To end the irrigation process, loosen the strap and slide it to the collection bag.

G – Oral stool modification

The successful use of the stool drainage system depends on the consistency of the stool to be drained being sufficiently able to flow. Stools that will flow can be controlled by means of oral stool modification prescribed for the patient. Any shifts in the patient's fluid balance should be monitored if the stools are consistently able to flow or liquid. Oral stool modification measures fundamentally require a medical prescription.

H – Sampling

Position the red closure strap approximately 5 cm below the sampling port and tighten the strap there. Insert the hub of the sampling syringe as far as possible into sampling port and collect the stool sample. Loosen the strap and slide it to the collection bag.

I – Changing and disposing of the collection bag

1. The collection bag should not remain connected to the draining catheter for more than 48 hours.
2. Avoid contamination with stools when changing the collection bag. Position the closure strap immediately adjacent to the connector and tighten it in such a way that it seals.
3. Seal the used bag using the associated closure cap. Dispose of the bag properly in accordance with hospital regulations. Comply with the specifications on the disposal of contaminated products/waste and for stools containing pathogens from patients with infectious diseases.
4. Connect the drainage system connector to the fresh stool collection bag. Twist the locking ring until it snaps into place.
5. Loosen the strap and slide it to the collection bag.

J – Removal and disposal of the drainage system

1. Completely deflate the drainage head and remove the drainage system carefully from the rectum.
2. Dispose of the used drainage system in accordance with local hygiene regulations. Comply with the specifications on the disposal of contaminated products/waste and for stools containing pathogens from patients with infectious diseases.

K – Storage conditions

Ambient temperature: 15 °C to 30 °C

There are no further storage-related restrictions for hygh-tec® odor-solution and hygh-tec® bag stool collection bags.

hygh-tec® odor-solution | product content

Stool drainage system including strap, stool collection bag, inflating syringe, instructions for use

The hygh-tec® bag stool collection bag is only used in conjunction with the hygh-tec® odor-solution stool drainage system.

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